



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
ABS Bedside Cabinet

Product Code :
SAMCOMEDICAL10021A003

Description :

Bedside Cabinet...ABS

Technical Specification :

Bedside Cabinet...ABS ABS Bedside Cabinet,attractive Blue color. Size: 475 ½ 470 ½ 755 mm. 1 drawer and 1 storage cabinet. Hidden pull-out tray for cups & articles. 2 towel holders (one on each side). Note: Available in Beige/Green color.



Samco Medical

Website: www.samcomedical.com,