



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Bedside Cabinet

Product Code :
SAMCOMEDICAL10021A006

Description :

Bedside Cabinet

Technical Specification :

Bedside Cabinet Size: 490 ½ 400 ½ 755 mm. With stainless steel top & bottom. Remaining framework made of epoxy coated mild steel. 1 drawer and 1 storage cabinet. 2 towel holders (one on each side).



Samco Medical

Website: www.samcomedical.com,