



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Bedside Locker std. Hospital Model

Product Code :
SAMCOMEDICAL10021A009

Description :

Bedside Locker std. Hospital Model

Technical Specification :

Bedside Locker std. Hospital Model Overall size : 40L x 40W x 82H cms. Machine pressed sheet box with door of heavy gauge. Assembled to CRC tubular legs of 2.5 cms. dia. fitted with PVC stumps. Stainless steel polished top of 22G. Pre-treated and epoxy powder coated.



Samco Medical

Website: www.samcomedical.com,