



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Examination Table Bed

Product Code :
SAMCOMEDICAL10023A002

Description :

Samcomedical Examination Table Bed

Technical Specification :

Examination Table /Bed Examination Table / Couch, 2 sections. Size: 1900 620 700 mm. Back section: 695 620 mm Back rest adjustable by double ratchet. Framework made of Stainless Steel. Mounted on protective Stumps. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,