



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Examination Table /Bed, General

Product Code :
SAMCOMEDICAL10023A003

Description :

Examination Table /Bed, General

Technical Specification :

Examination Table /Bed, General Examination Table / Couch, single section Size: 1900 i_l½ 620 i_l½ 700 mm. Framework made of Mild Steel Epoxy coated finish. Mounted on protective stumps. Freight saving knock-down construction.



Samco Medical

Website: www.samcomedical.com,