



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Examination Table - (With Accessories)

Product Code :
SAMCOMEDICAL10023A004

Description :

Examination Table - (With Accessories)

Technical Specification :

Examination Table - (With Accessories) Examination table, 2 sections, adjustable head section with ratchet mechanism to varying position. Legs mounted on rubber tips. Size : Table assembled : Approx. 1800L x 550W x 800H mm. Tubes: Approx. dia. 30 x 1.5 mm. Material : Resistant to corrosion (suitable for tropical environment) and decontamination (chlorine solution). Frame: Epoxy coated tubular steel. Mattress-pad : Foam density approx. 21 - 25 kg / m³. Fixed mattress-pad : Approx. 30H mm. Carrying capacity: Approx. 150 - 180 kg. Supports for leg holders with solid mechanism to fix leg holders at requested height and width. Knockdown construction



Samco Medical

Website: www.samcomedical.com,