



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Examination Tables

Product Code :
SAMCOMEDICAL10023A005

Description :

Examination Tables

Technical Specification :

Examination Table Examination table, 2 sections, adjustable head section with lever mechanism to varying position. Legs mounted on rubber tips. Size : Table assembled : Approx. 1800L x 550W x 800H mm. Framework made of mild steel epoxy coated finish with upholstered top. Carrying capacity : Approx. 150 - 180 kg. Knockdown construction



Samco Medical

Website: www.samcomedical.com,