



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Examination Table

Product Code :
SAMCOMEDICAL10023A006

Description :

Samcomedical Examination Table

Technical Specification :

Examination Table Examination table, 3 sections, with adjustable head and leg section with lever mechanism to varying position. Legs mounted on rubber tips. Size : Table assembled: Approx. 1800L x 550W x 800H mm. Framework made of chrome plated rectangular steel tubes with upholstered top. Carrying capacity: Approx. 150 - 180 kg. Knockdown construction.



Samco Medical

Website: www.samcomedical.com,