



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Stretcher Single Fold With Two Wheels

Product Code :

SAMCOMEDICAL10025A010

Description :

Samcomedical Stretcher Single Fold With Two Wheels

Technical Specification :

Stretcher Single Fold With Two Wheels Light Weight, easy clean, anodized aluminium alloy framework with 2 wheels. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 186x52x19cm. Folded Size (LxWxH): 93x52x10cm. Self-weight: 5 kg, Load Bearing: 159 kg



Samco Medical

Website: www.samcomedical.com,