



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Stretcher Single Fold With Adjustable Back Rest

**Product Code :**  
SAMCOMEDICAL10025A011

**Description :**

Stretcher Single Fold With Adjustable Back Rest

**Technical Specification :**

Stretcher Single Fold With Adjustable Back Rest Light Weight, easy clean, anodized aluminium alloy framework with adjustable back rest & 2 wheels. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 185x54x24cm. Folded Size (LxWxH): 92x54x14cm. Self-weight: 7 kg, Load Bearing: 159 kg



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),