



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Single Fold

Product Code :
SAMCOMEDICAL10025A012

Description :

Stretcher Single Fold

Technical Specification :

Stretcher Single Fold Light Weight, easy clean, anodized aluminium alloy framework. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 185x52x5cm. Folded Size (LxWxH): 93x52x6cm. Self-weight: 3 kg, Load Bearing: 159 kg



Samco Medical

Website: www.samcomedical.com,