



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Single Fold With Aluminium Top

Product Code :
SAMCOMEDICAL10025A013

Description :

Stretcher Single Fold With Aluminium Top

Technical Specification :

Stretcher Single Fold With Aluminium Top Light Weight, easy clean, anodized aluminium alloy framework. Washable, water proof & anti-fungus aluminium top with safety belts. Unfolded Size (LxWxH): 185x52x5cm. Folded Size (LxWxH): 94x50x7cm. Self-weight: 7 kg, Load Bearing: 159 kg



Samco Medical

Website: www.samcomedical.com,