



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Army 3 Fold With Two Wheels

Product Code :
SAMCOMEDICAL10025A014

Description :

Stretcher Army 3 Fold With Two Wheels

Technical Specification :

Stretcher Army 3 Fold With Two Wheels Light Weight, easy clean, anodized aluminium alloy framework with 2 wheels. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 188x45x22cm, Folded Size (LxWxH): 46x45x15cm. Self-weight: 8 kg, Load Bearing: 159 kg. Supplied with shoulder bag



Samco Medical

Website: www.samcomedical.com,