



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Double Fold

Product Code :
SAMCOMEDICAL10025A016

Description :

Stretcher Double Fold

Technical Specification :

Stretcher Double Fold Light Weight, easy clean, anodized aluminium alloy framework. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 215x54x15cm, Folded Size (LxWxH): 108x17x9cm. Self-weight: 5 kg, Load Bearing: 159 kg. Supplied with carry bag



Samco Medical

Website: www.samcomedical.com,