



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher 4 Fold

Product Code :
SAMCOMEDICAL10025A017

Description :

Stretcher 4 Fold

Technical Specification :

Stretcher 4 Fold Light Weight, easy clean, anodized aluminium alloy framework. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 215x55x14cm. Folded Size (LxWxH): 56x22x20cm. Self-weight: 8 kg, Load Bearing: 159 kg. Supplied with carry bag



Samco Medical

Website: www.samcomedical.com,