



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Stretcher 4 Fold

Product Code :
SAMCOMEDICAL10025A019

Description :

Samcomedical Stretcher 4 Fold

Technical Specification :

Stretcher 4 Fold Light Weight, easy clean, anodized aluminium alloy framework. Washable, water proof & anti-fungus fabric top. Unfolded Size (LxWxH): 210x53x12cm. Folded Size (LxWxH): 61x30x12cm. Self-weight: 6 kg, Load Bearing: 159 kg. Supplied with shoulder bag.



Samco Medical

Website: www.samcomedical.com,