



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Vacuum Mattress

Product Code :
SAMCOMEDICAL10025A020

Description :

Stretcher Vacuum Mattress

Technical Specification :

Stretcher Vacuum Mattress Made of high quality & high resistant PVC material. Adjust perfectly to the shape of body. It allows radiological X-Ray. Unfolded Size (LxWxH): 200x98x7cm, Folded Size (LxWxH): 30x30x50cm. First aid personnel can rapidly and safely adjust the mattress to be soft or hard by using double action pump. Self-weight: 5.5 kg, Load Bearing: 159 kg



Samco Medical

Website: www.samcomedical.com,