



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Stretcher Cum Wheel Chair Single Fold With 2 Wheels

**Product Code :**

SAMCOMEDICAL10025A021

**Description :**

Stretcher Cum Wheel Chair Single Fold With 2 Wheels

**Technical Specification :**

Stretcher Cum Wheel Chair Single Fold With 2 Wheels Designed for easily carrying patients in narrow and hard passing places. Light Weight, easy clean, anodized aluminium alloy framework. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 210x55x12cm, Folded Size (LxWxH): 102x55x16cm Folded in Wheelchair Back-H:138cm Seat-H:53cm. Self-weight: 9 kg, Load Bearing: 159 kg



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),