



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Stretcher Army Single Fold, with Telescopic Lifting Handles

Product Code :

SAMCOMEDICAL10025A022

Description :

Stretcher Army Single Fold, with Telescopic Lifting Handles

Technical Specification :

Stretcher Army Single Fold, with Telescopic Lifting Handles Light Weight, easy clean, aluminium alloy framework with telescopic lifting handles. Washable, water proof & anti-fungus fabric top with safety belts. Unfolded Size (LxWxH): 220x54x15cm. Folded Size (LxWxH): 180x12x15cm. Self-weight: 8 kg, Load Bearing: 159 kg.



Samco Medical

Website: www.samcomedical.com,