



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Camping Bed / Battlefield Stretcher

Product Code :
SAMCOMEDICAL10025A025

Description :

Camping Bed / Battlefield Stretcher

Technical Specification :

Camping Bed / Battlefield Stretcher Light Weight, easy clean, aluminium alloy framework. It is mainly used for battle field and outdoor carrying patients and wounded person. When unfolded, it can be used as a simple bed in the tent. Unfolded Size : 190L x 65 W x 42H cm Folded size : 93L x 20W x 20H cm Self Weight : 8 Kg. Load bearing : 159 Kg. Self Weight : 8 Kg



Samco Medical

Website: www.samcomedical.com,