



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Salvage Stretcher Yellow

Product Code :
SAMCOMEDICAL10025A026

Description :

Salvage Stretcher Yellow

Technical Specification :

Salvage Stretcher Yellow Designed for fractured and immobile serious patients. Light Weight, easy clean, anodized aluminium alloy framework with safety belts. Size (LxWxH): Max. unfolded Size (LxWxH): 214x42x7cm, Min. unfolded Size (LxWxH): 167x42x7cm. Detachable and folding. This Salvage stretcher can be used with head immobilizer. Length of stretcher adjustable according to patient's body. The foot-end is designed as a narrow-frame structure. Self-weight: 8 kg, Load Bearing: 159 kg



Samco Medical

Website: www.samcomedical.com,