



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Head Immobilizer

Product Code :
SAMCOMEDICAL10025A029

Description :

Head Immobilizer

Technical Specification :

Head Immobilizer Made of high-density plastic material, large ear holes for monitoring the patient's ear canal. It is used with Salvage stretcher & Spine Stretchers. Minimal interference with X-ray, MRI or CT scanning procedures. Waterproof plastic coating makes it easy to clean & prevents bacterial growth on it



Samco Medical

Website: www.samcomedical.com,