



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Blood Donor Chair Cum Bed

Product Code :

SAMCOMEDICAL10026A003

Description :

Samcomedical Blood Donor Chair Cum Bed

Technical Specification :

Blood Donor Chair Cum Bed Framework made of mild steel epoxy coated finish. In folding position, it can be used as chair. When chair is spread, it can be used as a bed and the patient can lie on it. Size: 1900 x 530 x 240/400 mm



Samco Medical

Website: www.samcomedical.com,