



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Baby Bassinet

Product Code :
SAMCOMEDICAL10020A002



Description :

Frame work made of ERW tubular pipe of 2.5 cms. dia. mounted on four 5 cms dia. castors. Removable CRIB of Prepex transparent Crib with mattress and R-trend positions. Trendlenburg / Rev. Trendlenburg. Laminated cupboard for storage, I.V Rod. Pre-treated and epoxy powder coated.



Samco Medical

Website: www.samcomedical.com,