



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Bowl Stand (Single)

Product Code :
SAMCOMEDICAL10017A005

Description :

Samcomedical Bowl Stand (Single)

Technical Specification :

Medical Bowl Stand (Single) ERW steel tubular pipe construction mounted on heavy polypropylene base with castors. Pre-treated and epoxy powder coated. Supplied complete with S.S. basin of 35 cm. dia. Knockdown construction.



Samco Medical

Website: www.samcomedical.com,