



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Bowl Stand (Double)

**Product Code :**  
SAMCOMEDICAL10017A007

**Description :**

Samcomedical Bowl Stand (Double)

**Technical Specification :**

Medical Bowl Stand (Double) ERW steel tubular pipe construction mounted on heavy polypropylene base with castors. Pre-treated and epoxy powder coated. Supplied complete with S.S. basin of 35 cm. dia.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),