



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Baby Cradle On Stand Swinging Type

Product Code :
SAMCOMEDICAL10020A001

Description :

MS tubular frame work mounted on strong rubber shoes. Stand Fitted with PVC stumps Wire mesh hanging cradle Pre-treated and Epoxy Powder coated. Size : 36" x 15" x 39" , Height upto CRIB top 98 cms. Optional Accessories: 50 mm dia castors.



Samco Medical

Website: www.samcomedical.com,