



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Scalp Vein Set

Product Code :
SAMCOMEDICAL500001A021



Description :

Scalp Vein Set Scalp Vein Set is used for rapid venous access and for the greater patient comfort during infusion. Manufactured from non toxic, medical grade PVC compound. Thin wall, short bevel Needle facilitates non traumatic penetration. Wings are for better handling and attachment to the skin. Wings are colour coded for instant size identification. Proximal end of the set is provided Female luer fitting for easy connection with the infusion set.



Samco Medical

Website: www.samcomedical.com,