



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Pediatric Bed With Side Railings

Product Code :
SAMCOMEDICAL10020A004



Description :

Pediatric Bed With Side Railings Made of Tubular Pipe framework. Top of M.S. Sheet with perforations. Built in back rest. Mounted on four 10 cms. dia castors (2 with brakes). Drop side, safety railing at both side. Pre-treated and epoxy powder coated finish. Size : 1370 L x 760 W cm x 610 H mm Pediatric Bed With Side Railings, Infant Bed & Cradle Exporters, Infant Bed & Cradle Manufacturers



Samco Medical

Website: www.samcomedical.com,