



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Blood Transfusion Set

Product Code :
SAMCOMEDICAL500001A025



Description :

Blood Transfusion Set Cylindrical drip chamber with and without vent provided to prevent passage of any clot into the patient. 18 G Needle for smooth flow. Moulded rubber with filter flash ball for extra medication. Efficient roller controller for accurate adjustment of infusion rate. E.T.O. Sterile and Pyrogen Free. Luer lock connectors also available



Samco Medical

Website: www.samcomedical.com,