



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Intravenous Transfusion Set

**Product Code :**

SAMCOMEDICAL500001A026

**Description :**

Intravenous Transfusion Set Extra strong, spike with/without air vent bacteria barrier at air inlet. Flexible drip chamber facilitate rapid adjustment of fluid level. Disc filter in drip chamber. Kink resistant non-toxic medical grade P.V.C. tubing. Strong roller clamp with a provision to hold the tube for changing the fluid container without disturbing the flow rate. Extra thick self sealing latex bulb. Luer lock connectors also available. E.T.O. Sterile and Pyrogen Free



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),