



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Examination Gloves Latex, Non Sterile, Disposable

Product Code :
SAMCOMEDICAL500001A029



Description :

Examination Gloves Latex, Non Sterile, Disposable Made from pure natural rubber Latex for greater finger dexterity. Ambidextrous. Non sterile. Disposable, for single use only. Slightly Powdered. Micro rough-finish for better grip. Weight 5 gram per pcs. Packing of 100 pcs per box.



Samco Medical

Website: www.samcomedical.com,