



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Pediatric Bed S.Steel

Product Code :
SAMCOMEDICAL10020A005



Description :

Semi-Fowler Bed with mechanically operated Back rest by crank mechanism. Steel frame work with 2 section top. S.S. side railings with Up & Down movement. With S.S. Shelf at the bottom for keeping blankets, basin etc. Tubular Head & Foot Bows also made of stainless steel. Mounted on four silent swivel castor (all with brakes). Provision for I.V. Rod. Freight saving knock-down construction. Size: 1840 \times 800 \times 580 mm. Pediatric Bed S.Steel, Infant Bed & Cradle Exporters, Infant Bed & Cradle Manufacturers



Samco Medical

Website: www.samcomedical.com,