



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Male External Catheter

**Product Code :**  
SAMCOMEDICAL500001A034



**Description :**

Male External Catheter Position Over Tip of Penis Then Unroll the latex sleeve over entire length. Use adhesive strap around top of sleeve to secure to penis avoiding constriction. Afix distal end to drainage tubing. Used for Urine incontinence in male patients. Made from natural latex rubber. Sterile, Individually packed in peelable pouch pack



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),