



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Urine Collection Bag (2000ml, Superior)

Product Code :
SAMCOMEDICAL500001A036



Description :

Urine Collection Bag (2000ml, Superior) Used for short & long term urine drainage. Non return valve to avoid the back flow of urine. Graduated with easy to read scale. Also available with top bottom outlet & moulded handle for easy handling and carrying facility. Sterile, Individually packed in Tear H.M pack.



Samco Medical

Website: www.samcomedical.com,