



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Urine Collection Bag with Measured Volume Meter

Product Code :
SAMCOMEDICAL500001A037



Description :

Urine Collection Bag with Measured Volume Meter With measured volume meter for measurement of urine input. 2000ml Urine Bag and 250ml Urine Bottle. Urine Bottle is directed to urine bag as to avoid any possibility of contamination. Push/Pull type bottom outlet for convenient emptying of urine bag. Sterile, Individually packed in peelable pouch pack.



Samco Medical

Website: www.samcomedical.com,