



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Nelaton Catheters

**Product Code :**  
SAMCOMEDICAL500001A041



**Description :**

Nelaton Catheters Used for bladder catheterisation through urethra. Made from non-toxic, non-irritant, medical grade PVC compound with radio opaque line. Distal end with 2 lateral eyes for maximum drainage & proximal end with funnel shape connector for easy connection to urine bag. Sterile, Individually packed in peelable blister pack.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),