



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Dressing Trolley Steel with Drawer

**Product Code :**  
SAMCOMEDICAL10015A001

**Description :**

Trolley made of stainless steel. Size: 640 x 420 x 800 mm. 2 S.S. shelves with guard rails on 3 sides. With 1 wastebin & 1 S.S. Bowl. S.S. Dressing tray attached at the top. 1 S.S. drawer below the first shelf. 4 swivel castor (2 with brakes).



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),