



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Endobronchial Suction Catheters

Product Code :
SAMCOMEDICAL500001A046



Description :

Endobronchial Suction Catheters Suction Catheters are used for removal of secretion from mouth, trachea and bronchial tube. Manufactured from non toxic, non irritant medical grade PVC Compound. Distal end is open with one lateral eye for effective drainage. Proximal end is provided with connector for easy connection to the suction equipment. Funnel connector are colour coded for instant size identification. Provided with vacuum control facility in finger tip and thumb control models. Sterile, individually packed in peelable blister pack



Samco Medical

Website: www.samcomedical.com,