



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Tracheostomy Tube

Product Code :
SAMCOMEDICAL500001A047



Description :

Tracheostomy Tube Tracheostomy tube is used to provide unobstructed airway to lungs by passing the larynx and upper respiratory tract. Manufactured from non toxic, non irritant medical grade PVC Compound readily conforming to body contours. Soft flexible wing fixed at proximal end helps in easy fixation. Terminal end is provided with 15mm connector. Sterile, individually packed in peelable pouch pack.



Samco Medical

Website: www.samcomedical.com,