



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Instrument Trolley

Product Code :
SAMCOMEDICAL10015A003

Description :

Instrument Trolley Frame work made of ERW tubular pipe mounted on four 7.5 cms. dia. castors. Two fixed stainless steel shelves of heavy gauge. Pre-treated and epoxy powder coated. Also available : Instrument Trolley - All stainless steel (NMTOY10002A-0086) Instrument Trolley Without railing. (NMTOY10002B-0087) Instrument Trolley Size: 68L x 45W x 82H cms. Instrument Trolley Size: 76L x 45W x 82H cms. Instrument Trolley Size: 91L x 45W x 82H cms.



Samco Medical

Website: www.samcomedical.com,