



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Endotracheal Tube

**Product Code :**  
SAMCOMEDICAL500001A053



**Description :**

Endotracheal Tube Used to provide a direct and unimpeded airway for gases to pass to and from the lungs. Manufactured from non toxic, non irritant medical grade PVC Compound. Tip is suitable for both Nasal and Oral intubation. Radio opaque line provided throughout the length for X-ray visualization. Sterile, individually packed in peelable pouch pack



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),