



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Umbilical Catheter

Product Code :
SAMCOMEDICAL500001A061



Description :

Umbilical Catheter Used to access the umbilical vessels of pre-matured infants or newly born babies. Manufactured from non toxic, medical grade PVC. Radio opaque line provided throughout the tube for checking the position of the catheter radiologically. Marking is provided at every cm to ascertain depth of insertion. Proximal end of catheter is provided with Luer Female connector for easy connection to the I.V. Therapy equipment. Sterile, individually packed in peelable pouch pack



Samco Medical

Website: www.samcomedical.com,