



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Scalpels, Disposable, Gamma Sterile

Product Code :
SAMCOMEDICAL500001A082



Description :

Scalpels, Disposable, Gamma Sterile Scalpels - Disposable, gamma sterile, stainless steel blade, factory fitted, in convenient graduated plastic handle (with scale). Blade protected by plastic cap. Individually in peel pack using high quality medical grade paper. Box of 10 pcs



Samco Medical

Website: www.samcomedical.com,