



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Salem Sump Tube / Gastro Duodenal Dual Flow Tube /  
Irrigation Tube

**Product Code :**  
SAMCOMEDICAL500001A110



**Description :**

Salem Sump Tube / Gastro Duodenal Dual Flow Tube / Irrigation Tube Dual flow type tube is used to prevent the accumulation of gastric secretions and air. Double lumen tube throughout the length provides air vent when the stomach has been evacuated of accumulated liquid, air or small particulate matter. Distal end is coned with ten lateral eyes. Radio opaque tube is marked at 45, 55, 65 and 75 m from the tip for accurate placement. Proximal end is provided with Universal Funnel Connector along with Double Taper connector for easy extension. Sterile, Individually packed in peelable pack



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),