



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Kehr T-Tubes

Product Code :
SAMCOMEDICAL500001A111



Description :

Kehr T-Tubes Kehr s T tube is used for short term post operative drainage of common bile duct. Manufactured from soft non toxic, medical grade PVC compound. Radio opaque line provided throughout the tube for X-ray visualization. Super smooth tube surface prevent clogging. Sterile, Individually packed in peelable pouch pack



Samco Medical

Website: www.samcomedical.com,