



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hospital Bedside Locker

Product Code :
SAMCOMEDICAL10021A001

Description :

Hospital Bedside Locker

Technical Specification :

Machine pressed double bent CRCA sheet. Super imposed stainless steel top with raised edge on three sides and one drawer MS tubular legs with rubber feet. Pre-treated and epoxy powder coated finish. Size(in) : 16 L x 16 W x 32 H Hospital Bedside Locker, Bedside Lockers Exporters, Bedside Lockers Manufacturers



Samco Medical

Website: www.samcomedical.com,