



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Ryles Tubes (Stomach Tube)

**Product Code :**  
SAMCOMEDICAL500001A112



**Description :**

Ryles Tubes (Stomach Tube) Ryle s tube is used for nutritional feeding purposes or aspiration of intestinal secretions. Manufactured from non toxic, medical grade PVC compound. Distal end is coned with steel balls sealed into the tube to facilitates easy insertion. Four lateral eyes for efficient drainage. The tube is marked at 50, 60 and 70cm from the distal end for accurate placement into the abdomen. Radio opaque line provided throughout the tube for X-ray visualization. Proximal end is provided with Universal Funnel Connector for easy extension. Sterile, Individually packed in peelable pouch pack



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),