



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Rectal Catheter

**Product Code :**  
SAMCOMEDICAL500001A114



**Description :**

Rectal Catheter Used for introduction of enema solution into rectum. Manufactured from non-toxic, non irritant medical grade PVC compound Distal end is open with one lateral eye for efficient drainage. Proximal end is provided with Funnel shape connector for easy connection to the suction equipment. Funnel connectors are colour coded for instant size identification. Sterile, Individually packed in peelable pack



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),