



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Infant Feeding Tubes, Sterile

**Product Code :**  
SAMCOMEDICAL500001A116



**Description :**

Infant Feeding Tubes, Sterile Infant Feeding tube is used for neonates and paediatric feeding. Manufactured from non toxic, medical grade PVC compound. Distal end is smooth coned shape for non traumatic intubation. Two lateral eyes for efficient drainage. The tube is marked at 20 cm from the distal end for accurate placement into the abdomen Radio opaque line provided throughout the tube for X-ray visualization. Proximal end is provided with luer mount for connection with feeding device. Sterile, Individually packed in peelable blister pack



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),